

**ISTITUTO ZOOPROFILATTICO SPERIMENTALE DELLE VENEZIE
EUROPEAN UNION REFERENCE LABORATORY FOR AVIAN INFLUENZA AND
NEWCASTLE DISEASE**



**SOP VIR 1004
HA AND NA SUBTYPING OF AVIAN INFLUENZA VIRUS BY
REAL-TIME RT-PCR
(Hoffmann et al., 2016; James et al., 2018)**

Standard operating procedure used by the EURL for AI and ND at the Istituto Zooprofilattico Sperimentale delle Venezie. Released on 08/04/21.

Contents

1. Purpose and field of application	3
2. References	3
3. Safety	3
4. Materials	3
4.1. Reagents	3
4.2. Equipment	4
5. Procedure	4
5.1. Controls	4
5.2. Preparation of samples and storage	4
5.3. Isolation of nucleic acids	4
5.4. Real-time RT-PCR	5
5.5. Data analysis	6
6. Interpretation of results	6
6.1. Reliability of controls	6
6.2. Diagnostic samples	7
7. Characteristics of the method	7

1. Purpose and field of application

This protocol describes the procedure to subtype the HA and NA of avian influenza virus (AIV) in biological samples by real-time reverse transcription polymerase chain reaction (RRT-PCR).

The protocol targets segment 4 (HA) and segment 6 (NA) regions conserved within different AIV subtypes using multiple oligonucleotides sets, based on the assays developed by Hoffmann et al. (2016) and James et al. (2018) and their subsequent updates.

The H5, H7 and H9 subtype identification is not included in this document. Please refer to specific RRT-PCR assays available at SOPs VIR 014, 143, 144, 1001.

The procedure can be applied to RNA purified from AIV isolates (normally allantoic fluid), tissue homogenates, stool, tracheal and cloacal swabs properly collected and preserved.

2. References

- B. Hoffmann, D. Hoffmann, D. Henritzi, M. Beer, T.C. Harder. Riems influenza A typing array (RITA): an RT-qPCR-based low density array for subtyping avian and mammalian influenza A viruses. *Sci Rep* 6:27211, 2016. doi: 10.1038/srep27211;
- J. James, M.J. Slomka, S.M. Reid, S.S. Thomas, S. Mahmood, A.M.P. Byrne, J. Cooper, C. Russel, B.C. Mollett, E. Agyeman-Dua, S. Essen, I.H. Brown, S.M. Brookes. Development and application of real-time PCR assays for specific detection of contemporary avian influenza virus subtypes N5, N6, N7, N8 and N9. *Avian Dis* 63(sp1):209-218, 2018. doi: 10.1637/11900-051518-Reg.1;
- Commission Decision 2006/437/EC approving a Diagnostic Manual for avian influenza as provided for in Council Directive 2005/94/EC;
- OIE - World Organization for Animal Health, Terrestrial Manual, Chapter 3.3.4. Avian influenza (Version adopted in May 2015);
- SOP VIR 1000 - Sample preparation and nucleic acids isolation for the detection and typing of Avian influenza virus and Avian Orthoavulavirus type 1 by molecular methods.

3. Safety

Individual laboratories are responsible for ensuring that all the procedures described in this document are conducted under high safety standards, including awareness on chemical and biological risks. According to the risk hazard, either BSL2 or BSL3 facilities must be used. Safety rules at individual laboratories must be agreed with the biosecurity and biosafety officer and acknowledged by all the staff members involved.

4. Materials

4.1. Reagents

For commercial products, refer to the manufacturer's instructions for storage conditions.

- Nuclease-free water;
- Oligonucleotides resuspension buffer (e.g. TE pH 8.0);

- AgPath-ID One-Step RT-PCR Reagents (Applied Biosystems);
- Subtype-specific oligonucleotides sets (see Annex I).

4.2. Equipment

- General molecular biology laboratory equipment and consumables;
- CFX96 Touch Deep Well Real-time PCR Detection System (Biorad) (or equivalent real-time PCR platform).

5. Procedure

5.1. Controls

To ensure test reliability, the controls listed in Table 1 must be included in each run.

Control Type	Definition
Negative process control (NPC)	Sample containing no target organism that is processed along with the samples starting from the nucleic acids isolation phase
Positive template control (amplification control, PTC)	Sample containing a known amount of the target RNA of the searched AIV subtype and stored at $\leq -70^{\circ}\text{C}$, that is processed along with the samples starting from the RRT-PCR phase
Negative template control (NTC)	Sample containing all PCR reagents but no target RNA, that is processed along with the samples starting from the RRT-PCR phase

Table 1. Controls to ensure test reliability

This procedure is applied as a downstream method to samples which have tested positive for AI by a different molecular diagnostic assay. In case this procedure is applied as a frontline method from a new sample aliquot, a positive process control (PPC) (i.e. sample containing the target organism that is processed along with the samples starting from the nucleic acids isolation phase) can be used to assess the reliability of the analytical process.

5.2. Preparation of samples and storage

For samples preparation (allantoic fluid, tissues, stool, tracheal and cloacal swabs), see SOP VIR 1000.

Clinical samples and viral isolates must be stored at refrigerated temperature ($2-8^{\circ}\text{C}$) until the completion of the analysis. For long term conservation, store samples at $\leq -70^{\circ}\text{C}$.

Purified nucleic acids can be refrigerated for a few hours prior RRT-PCR, otherwise they must be stored at $\leq -70^{\circ}\text{C}$.

5.3. Isolation of nucleic acids

RNA can be isolated either by manual or automatic methods.

The IZSVe validated the following nucleic acids isolation systems:

- NucleoSpin RNA kit (Macherey-Nagel);
- QIAamp Viral RNA Mini Kit (Qiagen);
- QIASymphony DSP Virus/Pathogen Midi kit (Qiagen) on the QIASymphony SP instrument (Qiagen);

- MagMAX Pathogen RNA/DNA Kit (Applied Biosystems) on the KingFisher Flex Magnetic Particle Processor (ThermoFisher Scientific).

5.4. Real-time RT-PCR

The preparation of the master mix for the RRT-PCR reaction has to be carried out in a clean dedicated area. The reagents should be kept refrigerated throughout the entire preparation phase. For the whole procedure, the use of filter tips and nuclease-free plastics is recommended.

Subtype-specific oligonucleotides pre-mix preparation

Subtype-specific oligonucleotides sets can be assembled in advance as pre-mixes using TE pH 8.0, and stored until use. Final concentrations of each oligonucleotide in the reaction mix are detailed in Annex I.

Reaction mix preparation

Prepare a master mix volume sufficient for the number of samples to be tested in each run. The volumes reported in Table 2 are per single reaction.

Component	Initial concentration	Final concentration	µl per reaction
Nuclease-free water	-	-	4.5
2X RT-PCR Buffer	2X	1X	12.5
Subtype-specific oligonucleotides pre-mix	-	-	2
25X RT-PCR Enzyme Mix	25X	1X	1
Master mix minus template			20
Template			5
Total reaction volume			25

Table 2. Volume and concentration of the reaction mix components

The master mix should be mixed thoroughly and 20 µl solution per sample must be pipetted in real-time PCR strip/plate.

Template (5 µl) must be added in a separate room. To minimize the risk of cross-contamination, it is recommended to add controls and samples in the following order: NTC, NPC, nucleic acids isolated from diagnostic samples, PTC. In case of NTC, 5 µl of nuclease-free water should be added.

Cycling conditions

Place the strip/plate in the real-time PCR apparatus. Set up the loading scheme as well as the thermal profile as per Table 3, and select a FAM compatible detection channel for fluorescence acquisition.

Step	Temperature and time	Repetitions
Reverse transcription	45°C for 10 min	1
RT inactivation/initial denaturation	95°C for 10 min	1
Denaturation	95°C for 15 sec	45
Annealing (*)	56°C for 20 sec	
Extension	72°C for 30 sec	

Table 3. Real-time RT-PCR thermal profile. (*) Fluorescence acquisition has to be done during the annealing phase

5.5. Data analysis

Upon completion of the amplification reaction, amplification plots must be critically assessed. The baseline and the threshold can be set either automatically or manually. The threshold should be placed above the background fluorescence noise, across the exponential phase of all the amplification curves (corresponding to the early linear phase of the logarithmic view). The use of a PTC with a known target concentration can act as a calibrator to enable standardization of data analysis and an approximate estimation of the viral load.

6. Interpretation of results

6.1. Reliability of controls

Test reliability is assured if the controls yield the expected results, as reported in Table 4. In case of invalid results, the causative reason must be investigated and proper actions have to be taken.

Control	Expected Result	Action in case of invalid control
PPC (if used)	Positive, i.e. increase in fluorescence from the FAM fluorophore yielding a sigmoidal (or logarithmic) amplification curve at the expected Ct value	Repeat the analysis from the nucleic acids extraction and check the PPC stock
NPC	Negative, i.e. absence of fluorescence increase from the FAM fluorophore, with no sigmoidal (or logarithmic) amplification curve	Repeat the analysis from the nucleic acids extraction and check nucleic acids isolation reagents
PTC	Positive, i.e. increase in fluorescence from the FAM fluorophore yielding a sigmoidal (or logarithmic) amplification curve at the expected Ct value	Repeat the analysis from the RRT-PCR and check the PTC stock
NTC	Negative, i.e. absence of fluorescence increase from the FAM fluorophore, with no sigmoidal (or logarithmic) amplification curve	Repeat the analysis from the RRT-PCR and check RRT-PCR reagents

Table 4. Assessment of test reliability. Ct = threshold cycle

6.2. Diagnostic samples

Criteria for data interpretation of diagnostic samples and subsequent actions are reported in Table 5.

Result	Interpretation and action
Increase in fluorescence from the FAM fluorophore yielding a sigmoidal (or logarithmic) amplification curve with Ct \leq 36	Positive. Proceed with downstream molecular and/or virological analyses, if required
Negative, i.e. absence of fluorescence increase from the FAM fluorophore, with no sigmoidal (or logarithmic) amplification curve	Negative. No action has to be taken
Weak increase in fluorescence from the FAM fluorophore yielding a sigmoidal (or logarithmic) amplification curve with Ct between 36 and 45	Dubious. Repeat the analysis from the nucleic acids extraction starting from a new aliquot of the sample together with the original extract. Alternatively, another confirmatory method (i.e. with alternative oligonucleotides sets) can be applied to the nucleic acids extracts. If PCR inhibition is suspected, dilute the purified nucleic acids 1:10 with nuclease-free water and repeat the analysis from the RRT-PCR and/or repeat the analysis from the nucleic acids extraction, further diluting the original biological sample (SOP VIR 1000). In case of recurrent dubious results, it is responsibility of the Head of the Laboratory to identify the actions to be taken

Table 5. Interpretation of real-time RT-PCR diagnostic results

7. Characteristics of the method

The IZSVe assessed the performances of the method employing AIV samples and clinical specimens available at the IZSVe repository, as well as materials, equipment and procedures as described above.

For any technical query, please contact eurl.ai.nd.secretariat@izsvenezie.it

Annex I. Oligonucleotides sets targeting all the HA and NA subtypes; [1] Hoffmann et al., (2016); [2] James et al. (2018)

Target	Oligonucleotide	Sequence 5'-3'	Reference	Concentration in the reaction mix
H1	IAV-H1-115F	ACA CAA TAT GTA TAG GYT AHC ATG C	[1]	1.6 µM
	IAV-H1-199R	GAG TGT GTY ACT GTY ACA TTC TT		1.6 µM
	IAV-H1-147FAM	FAM-TCD ACM GAC ACT GTW GAC ACA GTA CTN GA-BHQ1		0.4 µM
H2	IAV-H2-470-F	GAC ACA GCA YAC RAC AAC TGG	[1]	1.2 µM
	IAV-H2-601-R	GTG TTG TTG TAT GAT CYT TTR GCA A		1.2 µM
	IAV-H2-522-FAM	FAM-CCN TCA TTC TTC AGG AAC ATG GTY TGG-BHQ1		0.4 µM
H3	IAV-H3-1667-F	TGG ATT TCC TTT GCC ATA TCA TG	[1]	1.2 µM
	IAV-H3-1784-R	ATR CAC TCA AAT GCA AAT GTT GCA		1.2 µM
	IAV-H3-1753-FAM	FAM-CTA ATG TTG CCT YTY TGG CAG GCC CAC AT-BHQ1		0.4 µM
H4	IAV-H4-1586F	GAC TCA RGG ATA CAA RGA CAT	[1]	1.2 µM
	IAV-H4-1599F	AAG GAC ATC ATY CTY TGG ATT TC		1.2 µM
	IAV-H4-1686R	ACA AGC CCA CAA AAT RAA GGC		1.2 µM
	IAV-H4-1696R	TTC CRT TYT GAC AAG CCC ACA A		1.2 µM
	IAV-H4-1628FAM	FAM-TCC ATA TCA TGC TTY TTR CTC GTT GC-BHQ1		0.4 µM
H6	IAV-H6-1666F	CTT GGT GTG TAT CAA ATY CTT GC	[1]	1.6 µM
	IAV-H6-1776R	CAT TGA RCC ATT TGA RCA CAT CCA		1.6 µM
	IAV-H6-1693FAM-MGB	FAM-TAT AGT ACG GTA TCG AGC AGY CT-MGB		0.4 µM
H8	IAV-H8-1604F	TAC AAA ATT CTY AGC ATY TAC AGT AC	[1]	1.6 µM
	IAV-H8-1677R	ATT ARA CCT CCA GCA AYC AGG A		1.6 µM
	IAV-H8-1654FAM	FAM-TGC CAA GCA RAG ACT GGC CGC CA-BHQ1		0.4 µM
H10	IAV-H10-991F	GTT GCT TGC WAC MGG AAT GAG	[1]	1.2 µM
	IAV-H10-1041F	GCC TGT TTG GDG CRA TAG C		1.2 µM
	IAV-H10-1122R	TTT TGR TGT CKG AAR CCA TAC CA		1.2 µM
	IAV-H10-1173R	ATA GCT GCY TGA GTA CTY TTG TA		1.2 µM
	IAV-H10-1092FAM	FAM-ACC ATY CCT TCC CAT CCR TTY TCT A-BHQ1		0.5 µM
H11	IAV-H11-1510F	ARG TYA GGA ATG GAA CAT ATG AYC A	[1]	1.6 µM
	IAV-H11-1723R	CAA ATG GTA CAT CTA CAT GAY CCA		1.6 µM

Target	Oligonucleotide	Sequence 5'-3'	Reference	Concentration in the reaction mix
H11	IAV-H11-1626FAM	FAM-ATT TAC AGC TGC ATY GCA AGY AGT CT-BHQ1	[1]	1 µM
H12	IAV-H12-1607F	AGC ATC TAC AGC AGT GTY GC	[1]	1.6 µM
	IAV-H12-1707R	CAG AAA GTA CAA CGA AYA TTT CCA		1.6 µM
	IAV-H12-1674FAM	FAM-CCG AAA ATG AAA CCC CCA ATA ATC ATG A-BHQ1		0.4 µM
H13	IAV-H13-1241F	ATT GAC AAA ATG AAT GGR AAY TAT GAY TC	[1]	1.6 µM
	IAV-H13-1372R	AAG AAG YTT DGC ATT RTA TGA CCA		1.6 µM
	IAV-H13-1304FAM	FAM-ATA AAY ATG CTY GCA GAY AGR ATA GAT GAY GC-BHQ1		0.4 µM
H16	IAV-H16-1589F	GGG ATA AAR TTG AAR ACT GAR GA	[1]	1.6 µM
	IAV-H16-1708R	ACT GCT RCA TGC CCA CAK TAT		1.6 µM
	IAV-H16-1635FAM-MGB	FAM-TTT AYA GYT GCA TTG CAA GCA G-MGB		0.4 µM
N1	IAV-N1-3-F	AGR CCT TGY TTC TGG GTT GA	[1]	2 µM
	IAV-N1-3-R	ACC GTC TGG CCA AGA CCA		2 µM
	IAV-N1-3-FAM	FAM-ATY TGG ACY AGT GGG AGC AGC AT-BHQ1		0.4 µM
N2	IAV-N2-1367F	AGT CTG GTG GAC YTC AAA YAG	[1]	2.4 µM
	IAV-N2-1488R	AAT TGC GAA AGC TTA TAT AGV CAT		2.4 µM
	IAV-N2-1444.1FAM-MGB	FAM-CCA TCA GGC CAT GAG CCT-MGB		0.4 µM
N3	IAV-N3-1348F	AAY AGT ATA GTT ACT TTC TGY GG	[1]	1.6 µM
	IAV-N3-1422R	CCA ATG TTR GAA CCA TCH GG		1.6 µM
	IAV-N3-1373FAM	FAM-TAR ACA ATG AAC CTG GAT CGG GVA A-BHQ1		0.4 µM
N4	IAV-N4-1335F	GAC YAG TGG TAG TAG YAT YGC	[1]	1.6 µM
	IAV-N4-1345F	AGT AGY ATT GCR TTY TGT GGT GTT		1.6 µM
	IAV-N4-1437R	AAA TYA CTT GTC TAT GTC AAA DGG		1.6 µM
	IAV-N4-1387FAM	FAM-TGG TCR TGG CCY GAT GGC GCT CT-BHQ1		0.4 µM
N5	IAV-N5-1322F	AAG AGA GRA CWA GCA TTT GGA C	[1]	1.2 µM
	IAV-N5-1353F	CTC CAC KGT RTT TTG TGG TGT		1.2 µM
	IAV-N5-1421R	GGA AGA ATT GCK CCA TCA YC		1.2 µM
	IAV-N5-1426R	CAA AKG GAA GAA TTG CKC CAT CA		1.2 µM
	IAV-N5-1375FAM	FAM-TCM AGT GAG GTC CCA GGR TGG TC-BHQ1		0.4 µM
	Mod-N5-1353F	TTC CAC TGT GTT TTG TGG KGT		[2]
	Mod-N5-1426R	CGA ATG GTA GAA TTG CTC CAT CA	0.6 µM	

Target	Oligonucleotide	Sequence 5'-3'	Reference	Concentration in the reaction mix
N5	Mod-N5-1375FAM	FAM-TCY AGT GAG GTC CCA GGA TGG TC-BHQ1	[2]	0.4 µM
N6	IAV-N6-14F-New	TGA ARA TGA ATC CAA ATC ARA AGA TAA	FLI update of [1]	0.8 µM
	IAV-N6-100R	GCG ATT CCT ATT AGC AGG CTT AC		0.8 µM
	IAV-N6-43FAM-New	FAM-TGC ATT TCA GCA ACA GGA GTA ACA CTA TC-BHQ1		0.4 µM
	Mod-N6-10F	AGG GTG AAA ATG AAT CCA AAT CA	[2]	1.6 µM
	Mod-N6-14F	TGA AAA TGA ATC CAA ATC ARA AGR TAA		1.6 µM
	Mod-N6-97R	AAT TCC TAT YAG CAG RCT TAC YAC		1.6 µM
	Mod-N6-43FAM	FAM-TGC RTT TCA GCM ACA GGA RTR ACA CTA TC-BHQ1		0.4 µM
N7	IAV-N7-1305F	GTT GAA TTA ATW AGA GGA AGR CC	[1]	1.6 µM
	IAV-N7-1430R	GAT YTG TGC CCC ATC RGG GA		1.6 µM
	IAV-N7-1383FAM	FAM-AGC CCA DTC YCA GTT GGG TCY GGT TC-BHQ1		0.4 µM
	Mod-N7-1305F	GTT GAA TTR ATT AGA GGR AGR CC	[2]	1.6 µM
	IAV-N7-1430R	GAT YTG TGC CCC ATC RGG GA		1.6 µM
	IAV-N7-1383FAM	FAM-AGC CCA DTC YCA GTT GGG TCY GGT TC-BHQ1		0.4 µM
N8	IAV-N8-1296F	TCC ATG YTT TTG GGT TGA RAT GAT	[1]	1.6 µM
	IAV-N8-1423R	GCT CCA TCR TGC CAY GAC CA		1.6 µM
	IAV-N8-1354FAM	FAM-TCH AGY AGC TCC ATT GTR ATG TGT GGA GT-BHQ1		0.8 µM
	Mod-N8-1296F	YCC CTG YTT TTG GGT CGA AAT GAT	[2]	0.8 µM
	Mod-N8-1423R	GCT CCA TCG TGC CAT GAC CA		0.8 µM
	Mod-N8-1354FAM	FAM-TCT AGT AGC TCC ATT GTA ATG TGT GGA GT-BHQ1		0.4 µM
N9	IAV-N9-1363F	AGY ATA GTA TCR ATG TGT TCC AG	[1]	1.6 µM
	IAV-N9-1439R	AAG TAC TCT ATT TTA GCC CCA TC		1.6 µM
	IAV-N9-1393FAM	FAM-TTC CTB GGA CAA TGG AAC TGG CC-BHQ1		0.4 µM
	IAV-N9-1363F	AGY ATA GTA TCR ATG TGT TCC AG	[2]	1.6 µM
	Mod-N9-1439R	AAG TAC TCT ATT YTA GCC CCR TC		1.6 µM
	IAV-N9-1393FAM	FAM-TTC CTB GGA CAA TGG AAC TGG CC-BHQ1		0.4 µM